



VIKTOR
SURGICAL INSTRUMENTS


2019

Dental, Veterinary and Manicure Instruments

SECTION
ENDOSCOPIC INSTRUMENTS


Grasping Forceps

465



Grasping Forceps

Straight, serrated jaws with proximal recess, double action.



465-3956
Ratchet, insulated shaft,
5mm, 32cm




465-3957
Monopolar,
5mm, 32cm

465-3958
Monopolar,
5mm, 45cm




465-3959
Roto-Lok, ratchet,
5mm, 32cm



465-3960
Roto-Lok, locking lever,
5mm, 32cm



465-3961
Roto-Cam, monopolar,
5mm, 32cm




Tc Grasping Forceps


Double action
Carb-Bite tips, 22mm long
jaws maintain firm grasp on
delicate tissue, double action.



465-3962
Roto-Lok, ratchet,
5mm, 32cm



465-3963
Roto-Lok, locking lever,
5mm, 32cm




Strong Tc Grasping Forceps

Heavy Carb-Bite inserts
allow jaws to maintain firm
hold on tissue, double action.




465-3964
Roto-Lok, ratchet,
5mm, 37cm



465-3965
Roto-Lok, locking lever,
5mm, 37cm


466

Grasping Forceps




Grasping Forceps

Straight, serrated jaws with proximal recess, double action.




465-3956
Ratchet, insulated shaft,
5mm, 32cm



465-3957 465-3958
Monopolar, Monopolar,
5mm, 32cm 5mm, 45cm




465-3959
Roto-Lok, ratchet,
5mm, 32cm



465-3960
Roto-Lok, locking lever,
5mm, 32cm



465-3961
Roto-Cam, monopolar,
5mm, 32cm




Tc Grasping Forceps


Double action
Carb-Bite tips, 22mm long jaws maintain firm grasp on delicate tissue, double action.



465-3962
Roto-Lok, ratchet,
5mm, 32cm



465-3963
Roto-Lok, locking lever,
5mm, 32cm




Strong Tc Grasping Forceps

Heavy Carb-Bite inserts allow jaws to maintain firm hold on tissue, double action.




465-3964
Roto-Lok, ratchet,
5mm, 37cm



465-3965
Roto-Lok, locking lever,
5mm, 37cm



Grasping Forceps

467




Retraction Grasping Forceps

Long-ridged jaws with 2 x 3 blunt interdigitating teeth for strong grasp of heavy tissue.



Atrau-Allis Grasping Forceps

Unique jaw design offers strong atraumatic grasp, double action



467-3974
Ratchet, insulated shaft,
5mm, 32mm

467-3975
Ratchet, insulated shaft,
5mm, 45mm




Atrau-Allis LJ Grasping Forceps


Unique jaw design in 20mm length offers strong atraumatic grasp, double action



467-3976
Roto-Lok, ratchet,
5mm, 37mm



467-3977
Roto-Lok, locking lever,
5mm, 37mm




467-3978
Double action,
Roto-Lok, locking lever,
5mm, 37mm



467-3979
Double action,
Roto-Lok, ratchet,
5mm, 37mm



467-3980
Roto-Lok, ratchet,
5mm, 32mm




467-3981
Roto-Lok, locking lever,
5mm, 32mm

467-3982
Roto-Lok, ratchet,
5mm, 37mm


467-3983
Roto-Lok, locking lever,
5mm, 37mm

468


Grasping Forceps



Allis Grasping Forceps
Interlocking distal teeth with fine-serrated inner jaws, double action.



468-3984
Ratchet, insulated shaft,
5mm, 32mm



Retrograde Grasping Forceps
Retrograde, deep-serrated teeth allow for non-slip hold, double action.



468-3985
Roto-Lok, ratchet,
5mm, 32mm



Grasping and Dissecting Forceps
Cup-shaped, toothed jaws, double action.



468-3986
Ratchet,
5mm, 32mm



Grasping and Extracting Forceps
Retrograde, deep-serrated teeth



468-3987
Monopolar
5mm, 32mm



468-3988
Spring-loaded
10mm, 32mm

Grasping Forceps

469



Claw Forceps

2x3 interlocking teeth for penetrating hold and extraction of tissue.



469-3989
Spring-loaded
10mm, 32mm




Claw Forceps

2x3 interlocking teeth for penetrating hold and extraction of tissue.



469-3990
Spring-loaded
5mm, 32mm




Debakey Forceps

Fenestrated 45mm long jaws with Debakey-type teeth, double action.



469-3991
Ratchet,
5mm, 32mm



Debakey Forceps


Fenestrated 55mm long jaws with Debakey-type teeth, double action.



469-3992
Roto-Lok, ratchet,
10mm, 37mm


470

Grasping Forceps




Paddle Babcock Forceps

Flat paddles with fine, straight Serrations for atraumatic hold of delicate tissue, double action.




470-3993
Roto-Lok, ratchet,
10mm, 37mm




Paddle Babcock Forceps

Flat paddles with fine, straight Serrations for atraumatic hold of delicate tissue, double action.



470-3994
Roto-Lok, ratchet,
5mm, 37mm




470-3995
Roto-Lok, locking lever,
5mm, 37mm




Babcock Forceps

Traditional style Babcock jaws designed with six rows of fine serrations for atraumatic, yet secure grasp of bowel and other delicate tissue, double action.



470-3996
Roto-Lok, ratchet,
10mm, 37mm




Babcock Forceps

Rounded jaws with 2x3 rows of De-Bakey-type teeth for atraumatic, yet secure grasp of bowel and other delicate tissue, double action



470-3997
Ratchet,
10mm, 32mm



470-3998
Roto-Lok, ratchet,
10mm, 32mm

Grasping Forceps

471



Babcock Forceps


Traditional design for grasping delicate tissue, double action.



471-3999
Ratchet
5mm, 32cm



471-4000
Roto-Lok, ratchet
5mm, 32cm




471-4001
Roto-Lok, locking lever,
5mm, 37cm




Swanstron Grasping Forceps


Unique design of 1x2 De-Bakey-type teeth with triangular fenestration provides an effective hold with atraumatic grasp, double action.



471-4002
Roto-Lok, ratchet,
5mm, 37cm



471-4003
Roto-Lok, locking lever,
5mm, 37cm




Duval Grasping Forceps

Triangular jaws with deep-serrated surfaces for non-slip, firm grasp of tissue, double action.



471-4004
Ratchet,
10mm, 32cm



471-4005
Roto-Lok ratchet,
10mm, 32cm

472

Grasping Forceps



Russian Grasping Forceps
40mm long jaws with heavy serration for effective tissue, double action.



472-4006
Ratchet,
10mm, 32cm



Foerster Sponge Forceps
45mm jaws firmly hold tissue, double action.




472-4007
Ratchet,
10mm, 32cm



Hunter² Clamps
Right-angle shaft with Atrau-Allis teeth for secure grasp and manipulation of tissue around structures, e.g stomach/esophagus in fundoplication, double action.



472-4008
Ratchet,
5mm, 32cm




Hunter² Bowel Grasper
Fenestrated, 20mm jaws for atraumatic manipulation of delicate tissue, e.g bowel stomach, double action.



472-4009
Roto-Lok,
5mm, 32cm


Grasping Forceps

473



Stone Extractor

To facilitate stone removal once neck of gallbladder has been externalized.




473-4010
19.0 cm / 7½"



Stone Extractor

Fenestrated, serrated, oval-shaped jaws that converge at distal end to prevent slippage, double action.



473-4011
Monopolar
5mm, 32cm

473-4012
Monopolar
5mm, 45cm



Strong Atraumatic Grasping Forceps

Fenestrated, serrated, curved jaws with split-steel construction, for fixation of ovary of manipulation of pelvic tissues.



473-4013
Spring handle
5mm, 32cm

473-4014
Spring handle
5mm, 45cm

474

Dissecting Forceps




Micro Dissecting Forceps

Delicate, fenestrated, serrated jaws for gentle blunt dissection, double action.



474-4015
Monopolar
5mm, 32cm



474-4016
Roto-Cam, monopolar
5mm, 32cm



Needle Nose Forceps

Slender, tapered and fine-serrated jaws.



474-4017
Monopolar
5mm, 32cm




Tapered Dissecting Forceps

Tips of serrated, tapered jaws meet at fine point to allow delicate, yet firm, grasp and dissection of thin adventitial tissue, double action.



474-4018
Monopolar
5mm, 32cm



474-4019
Roto-Cam, monopolar
5mm, 32cm



Micro Grasping Forceps

3mm small, tapered jaws for dissection of small structures.



474-4020
Monopolar
5mm, 32cm

Dissecting Forceps

475




Bullet-Nose Dissecting Forceps

Bullet-shaped, tapered, serrated jaws for dissection, double action.



475-4021
Roto-Cam, monopolar,
5mm, 32mm



Tapered Maryland Dissecting

Curved, finely-tapered, serrated jaws with wide spread, double action.



475-4022
Roto-Cam, monopolar,
5mm, 32mm




Maryland Dissecting

Curved, tapered, serrated jaws with wide spread, double action.



475-4023
Monopolar,
5mm, 32mm



475-4024
Roto-Cam, monopolar,
5mm, 32mm



Maryland Dissector

Slender, full-curved, long serrated jaws with wide spread, double action.



475-4025
Roto-Cam, monopolar,
10mm, 32mm

476

Dissecting Forceps and Hook



Mixer Dissecting Forceps

Slender, full-curved, long serrated jaws with wide spread, double action.



476-4026
10mm, 32mm




Mixer Spreader-Dissectors

Jaws design offers a true right-angle for dissection around delicate structures; double action. The distinctive ergonomic handle minimize hand fatigue while allowing the surgeon fine control for precise dissection.



476-4027
Roto-Cam, monopolar
10mm, 28cm



476-4028
Roto-Cam, monopolar,
right-handle, 10mm, 32cm




476-4029
10mm, 28cm




Mixer Spreader-Dissectors

Jaws design offers a true right-angle for dissection around delicate structures; double action. The distinctive ergonomic handle minimize hand fatigue while allowing the surgeon fine control for precise dissection.



476-4030
Roto-Cam, monopolar,
5mm, 28cm



476-4031
Roto-Cam, monopolar,
ring-handle, 5mm, 32cm



Appel³ Spreader-Dissecting Hook


Right-angle tip allows rapid blunt dissection and isolation of structures.




476-4032
10mm, 40cm

Dissecting Forceps


477



Supercut Scissors
 Razor-sharp upper blade with finely-serrated lower blade allows clean, precise cut without tissue slippage.



477-4033 5mm, 32cm 477-4034 5mm, 45cm



Supercut Endo Scissors
 Electrosurgical instrument with curved, Supercut blades and precise cutting action offers mechanical cutting, electrosurgical cutting, blunt dissection, electrosurgical coagulation- all in one reusable, cost-effective instrument; double action.



477-4035 Monopolar 5mm, 32cm 477-4036 Monopolar 5mm, 45cm



Hook Scissors
 Blades designed to lift, isolate and transect tissue.




477-4037 Roto-Cam, monopolar 5mm, 32cm



477-4038 Monopolar, 5mm, 32cm 477-4039 Monopolar, 5mm, 45cm 477-4040 Non-insulated, 5mm, 32cm



477-4041 Roto-Cam, Monopolar, 5mm, 32cm




Micro Scissors
 Fine, delicate, 3mm blades for small incision.



477-4042 Curved Left, 5mm, 32cm



477-4043 Straight, 5mm, 32cm



Mayo Operating Scissors
 Large, sturdy blades for incising heavy tissue.



477-4044 10mm, 32cm

478


Detach Scissors System



Detach Scissors Insert


478-4044

Cleaning brushes, pkg/6




Supercut Detach Scissors

Unique supercut blade has a finely-serrated lower blade that helps prevent tissue slippage, curved blades meet the contour of structures while the blunt-shaped distal tip allows tissue manipulation and cauterization with blades closed, double action.



Hook Detach Scissors

Classic hook blades designed to lift, isolate and transect tissue.



**Micro Detach Scissors
Curved**

Delicate, sharp-sharp, curved blades.

Detach Handle/Shaft

478-4045

Roto-Cam, monopolar,
5mm, 32cm



478-4046

Roto-Cam, monopolar,
5mm, 32cm

478-4047

Insert 5mm, 32cm

478-4048

Insert 5mm, 32cm, pkg/3



478-4049

Roto-Cam, monopolar,
5mm, 32cm

478-4050

Insert 5mm, 32cm

478-4051

Insert 5mm, 32cm, pkg/3



478-4052

Roto-Cam, monopolar,
5mm, 32cm

478-4053


Insert 5mm, 32cm

478-4054

Insert 5mm, 32cm, pkg/3

Dissecting Forceps


479



Micro Detach Scissors

Straight


Delicate, sharp-sharp, straight blades



479-4055
Roto-Cam, monopolar,
5mm, 32mm


479-4056
Insert, 5mm, 32mm

479-4057
Insert, 5mm, 32mm, plg/3



Metzenbaum Detach Scissors

Long, curved, tapered blades with wide spread, double action



479-4058
Roto-Cam, monopolar,
5mm, 32mm


479-4059
Insert, 5mm, 32mm

479-4060
Insert, 5mm, 32mm, plg/3



"Baby" Metzenbaum Detach Scissors


Shorter, curved, very tapered blades reach more narrow and delicate dissection areas, double action



479-4061
Roto-Cam, monopolar,
5mm, 32mm

479-4062
Insert, 5mm, 32mm

479-4063
Insert, 5mm, 32mm, plg/3



Mini Detach Scissors

Curved, tapered, delicate blades, double action



479-4064
Roto-Cam, monopolar,
5mm, 32mm

479-4065
Insert, 5mm, 32mm

479-4066
Insert, 5mm, 32mm, plg/3


480

Biopsy instruments



Biopsy Forceps

Two sharp teeth, one on each jaw, for biopsy or secure hold of thick-walled tissue.




480-4067
1x1 teeth,
spring-loaded,
5mm, 32mm

480-4068
1x1 teeth,
spring-loaded,
5mm, 45mm



Biopsy Punches

Serrated lower jaw.




480-4069
Monopolar,
5mm, 32mm

480-4070
Monopolar,
5mm, 45mm



Hook Biopsy Punches

Frangenheim-type electrosurgical biopsy punch with hook blade.



480-4071
Monopolar,
5mm, 32mm

480-4072
Monopolar,
5mm, 45mm



Biopsy Punches

Single tooth in upper jaw.




480-4073
Monopolar,
5mm, 32mm

480-4074
Monopolar,
5mm, 45mm



Biopsy Forceps

Oval-shaped fenestrated jaws.




480-4075
Monopolar,
5mm, 32mm

480-4076
Monopolar,
5mm, 45mm

Specialization

481




Micro Biopsy Forceps

3mm jaws for biopsy of small, delicate tissues.



481-4077
Spring-loaded,
5mm, 32cm



Spoon Forceps

Wide cup-shaped jaws for tissue or stone removal.



481-4078
Spring-loaded,
5mm, 32cm



Cholangiogram Clamp

Channel for placement of 5Fr or smaller cholangiogram catheter.



481-4079
Ratchet,
5mm, 32cm




Cholangiogram Guide

Curved tip for simple placement of cholangiogram catheter.



481-4080
5mm, 32cm



Bessler/Treat' Suture Forceps

Uniquely designed jaws allow easier grasp and positioning tubular structure for placement of purse string with straight of curved needle.



481-4081
5mm, 32cm



VIKTOR

SURGICAL INSTRUMENTS

Head Office:

533 Keswick Ct Rocheter Hills Michigan USA 48306

UAE Office:

NexGen Gadgets FZE, SM-Oce-C1- 518G, Ajman Free Zone

South East Asia Office:

324/34 Samprawuth Rd, Bangna, Bangna Bangkok 10260

Tel. : +662 107 4183

Mr.Tony +661 599 6754

E-mail : tony@viktormedical.com